# Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)

**Catalog No: PA0027**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## **Instructions for Use**

Please read before using this product.

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### **Check the integrity of the packaging before use.**

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# Bond™ Ready-To-Use Primary Antibody Melanoma Marker (HMB45)

Catalog No: PA0027

## Intended Use

*This reagent is for in vitro diagnostic use.*

Melanoma Marker (HMB45) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human HMB45 antigen in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using an automated Bond system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

## Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). Melanoma Marker (HMB45) primary antibody is a ready-to-use product that is recommended for use with either the Bond Polymer Refine Detection or Bond Polymer Refine Red Detection. The demonstration of human HMB45 antigen is achieved by first allowing the binding of Melanoma Marker (HMB45) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with an automated Bond system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

## Reagents Provided

Melanoma Marker (HMB45) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

## Clone

HMB45.

## Immunogen

Extract of pigmented melanoma metastases from lymph nodes.

## Specificity

Human HMB45 antigen.

## Subclass

IgG1, Kappa.

## Total Protein Concentration

Approx 10 mg/mL.

## Antibody Concentration

Greater than or equal to 1.3 mg/L as determined by ELISA.

## Dilution and Mixing

Melanoma Marker (HMB45) primary antibody is optimally diluted for use on the automated Bond system. Reconstitution, mixing, dilution or titration of this reagent is not required.

## Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the automated Bond system.

## Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Melanoma Marker (HMB45) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

### Instructions for Use

Melanoma Marker (HMB45) primary antibody is recommended for use on an automated Bond system in combination with either Bond Polymer Refine Detection (DS9800) or Bond Polymer Refine Red Detection (DS9390). The recommended staining protocols for Melanoma Marker (HMB45) primary antibody are IHC Protocol F when using Bond Polymer Refine Detection and IHC Protocol J when using Bond Polymer Refine Red Detection. Enzyme pretreatment is recommended using Bond Enzyme 1 for 5 minutes.

### Results Expected

#### Normal Tissues

Melanoma Marker (HMB45) showed positive granular, cytoplasmic staining in melanocytes in the basal layer of epidermis. All other normal tissues were negative. (Total number of cases stained = 99).

#### Tumor Tissues

Melanoma Marker (HMB45) stained 13/13 malignant melanomas. No staining was observed in a variety of other tumors. (Total number of cases stained = 101).

**Melanoma Marker (HMB45) is recommended for use as part of a panel of antibodies for the diagnosis of melanomas.**

### Product Specific Limitations

Melanoma Marker (HMB45) is recommended for use with either Bond Polymer Refine Detection or Bond Polymer Refine Red Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

### Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

### Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

### Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kapur RP, Bigler SA, Skelly M, et al. Anti-melanoma monoclonal antibody HMB45 identifies an oncofetal glycoconjugate associated with immature Melanoma Markers. The Journal of Histochemistry and Cytochemistry. 1992; 40(2):207-212.
5. Gown AM, Vogel AM, Hoak D, et al. Monoclonal antibodies specific for melanocytic tumours distinguish subpopulations of melanocytes. American Journal of Pathology. 1986;123(2):195-203.

ProClin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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ProClin™ 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.