ONETOUCH®

Ultra®

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ONETOUCH[®]
Ultra[®]

OneTouch[®] changes everything.[®]
OneTouch.[®] Ça change tout.[®]

OneTouch[®] Customer Care Line / Ligne InfoSoins OneTouch[®] :

Canada 1 800 663-5521

U.S.A./É.-U. 1 800 227-8862

Monday–Friday 9am–8pm Eastern Time 6am–5pm Pacific time

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Download from www.Somanuals.com. All Manuals Search on www.Somanuals.com.

ONETOUCH[®] *Ultra[®]*

DIABETES BLOOD GLUCOSE MONITORING SYSTEM

System Owner's Booklet

Dear OneTouch® Ultra® System Owner:

You have chosen one of the best blood glucose monitoring systems available. This booklet has important information you must know about the OneTouch® Ultra® System. Please read it carefully.

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that keeping blood glucose levels close to normal can reduce the risk of diabetes complications by up to 60%.* The results you get with OneTouch® Ultra® System can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

A warranty registration card is included with your OneTouch® Ultra® System. Please complete it and mail it to us; if you prefer to fill it out online, please visit www.OneTouch.ca

*American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Table of Contents

ABOUT YOUR NEW SYSTEMiv	USING THE METER MEMORY34
The Complete OneTouch® Ultra® Blood Glucose Monitoring System5	DOWNLOADING TEST RESULTS TO A PERSONAL COMPUTER36
OneTouch® Ultra® Blood Glucose Meter6	COMPARING METER AND LABORATORY RESULTS38
OneTouch® Ultra® Test Strips10	SETTING THE METER40
BEFORE TESTING12	CARING FOR YOUR OneTouch® Ultra® SYSTEM44
Checking the Display12	Meter44
Coding the Meter12	OneTouch® UltraSoft™ Sampler44
Checking the System15	Battery45
TESTING21	DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE48
Preparing the OneTouch® UltraSoft™ Adjustable Blood Sampler22	SPECIFICATIONS58
Fingertip Blood Sampling25	GUARANTEE59
Forearm Blood Sampling26	INDEX60
Step-by-Step Test Procedure30	
Used Lancet and Test Strip Disposal32	
Special Messages33	

About Your New System

The OneTouch® Ultra® System uses the latest blood glucose monitoring technology. It measures the glucose content of a blood sample by means of an electrical current produced in the test strip and sent to the meter for measurement.

Test results are “plasma-calibrated.” This makes it easier for you and your diabetes-care team to compare your meter results with laboratory tests.

Your OneTouch® Ultra® Blood Glucose Monitoring System consists of three main products: the OneTouch® Ultra® Blood Glucose Meter, OneTouch® Ultra® Test Strips (sold separately), and OneTouch® Ultra® Control Solution. These products have been designed, tested, and proven to work together to produce accurate blood glucose results. Use no other test strips or control solution with your meter.

The OneTouch® Ultra® System is intended for use outside the body (in vitro diagnostic use). It should be used only for testing fresh capillary whole blood samples for glucose (sugar). It should not be used for the diagnosis of diabetes or for testing newborns.

CAUTION: *Before using any product to test your blood glucose, read all instructions and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional. These recommendations apply to all blood glucose monitoring systems and are supported by the Diabetes Educator Section of the Canadian Diabetes Association.*

Important Information

- Severe dehydration resulting from excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- Test results below 3.3 mmol/L (60 mg/dL) mean low blood glucose (hypoglycemia). Test results greater than 10.0 mmol/L (180 mg/dL) mean high blood glucose (hyperglycemia). If you get results below 3.3 mmol/L (60 mg/dL) or above 10.0 mmol/L (180 mg/dL), and do not have symptoms, first repeat the test. If you have symptoms or continue to get results that fall below 3.3 mmol/L (60 mg/dL) or above 10.0 mmol/L (180 mg/dL), follow the treatment advice of your healthcare professional.
- If you are experiencing symptoms that are not consistent with your blood glucose test results AND you have followed all instructions described in the OneTouch® Ultra® Owner's Booklet, call your healthcare professional.
- A red blood cell count (hematocrit) that is either very high (above 55%) or very low (below 30%) can cause false results.

The Complete OneTouch® Ultra® Blood Glucose Monitoring System

- OneTouch® Ultra® Meter
- Owner's Booklet
- Quick Reference Guide
- OneTouch® UltraSoft™ Adjustable Blood Sampler
- Optional OneTouch® UltraClear™ Cap
- OneTouch® UltraSoft™ Sterile Lancets
- OneTouch® Ultra® Control Solution
- OneTouch® Ultra® Test Strips (Sold Separately)
- Carrying Case
- Warranty Registration Card
- Logbook
- One 3.0 V Lithium Battery (Installed)

OneTouch® Ultra® Blood Glucose Meter

DISPLAY

Symbols, simple messages, and test results appear here.

M BUTTON

Used to turn meter on to enter:

- ▶ setting mode
- ▶ memory mode

DATA PORT

Used to download your test results to a computer.

TEST PORT

Insert the OneTouch® Ultra®
Test Strip here.

C BUTTON

Used to:

- ▶ change date, time and code number
- ▶ indicate control solution tests
- ▶ review test results in memory

OneTouch® Ultra® Meter Display

CTL

Indicates a control solution test result.

CODE

Appears with the code number of the test strips.

BLOOD DROP SYMBOL

Tells you when to apply the sample.

DECIMAL POINT

Appears in test results only when unit of measure is set to mmol/L.

BATTERY SYMBOL


Warns when the battery is low or must be replaced.

KETONES?

Appears when a test result is above 13.3 mmol/L (240 mg/dL) to suggest ketone testing.

MEM

Indicates a test result stored in memory.

TEST RESULT AREA

Test results are displayed here.

MMOL/L


Unit of measure. Millimoles per litre (mmol/L) is the standard unit in Canada.

MG/DL

Unit of measure used in other parts of the world.


OneTouch® Ultra® Test Strips (Sold Separately)

The OneTouch® Ultra® System measures the amount of glucose in whole blood. Blood is applied to the TOP EDGE of the OneTouch® Ultra® Test Strip and is automatically drawn into the reaction cell where the reaction takes place.

Top Edge

Apply a drop of blood to the narrow channel here in the top edge of the test strip.

Confirmation Window

Check here to confirm if enough blood has been applied.

Contact Bars

Insert this end of the test strip, contact bars facing up, into the meter. Push it all the way in until it will go no further.

Important Test Strip Information


- Store test strip vials in a cool, dry place below 30°C (86°F). Keep away from direct sunlight and heat. Do not refrigerate.
- Store your test strips in their **original vial only**. To avoid damage or contamination, do not transfer test strips to any other place.
- Do not use test strips beyond the expiration date printed on the package since they may cause inaccurate results.
- After removing a test strip from the vial, replace the vial cap immediately and close it tightly.
- With clean, dry hands, you may touch the test strip anywhere on its surface.
- Use each test strip immediately after removing it from the vial.
- Count three months from the date you first open a new vial of test strips and write this date on the vial label. Throw test strips and vial away after this discard date.
- Apply only OneTouch® Ultra® Control Solution or a blood sample to the test strip.
- Do not bend, cut, or modify test strips in any way.
- OneTouch® Ultra® Test Strips are for single use only. **Never reuse a test strip that has had either blood or control solution applied to it.**
- Refer to additional information in the OneTouch® Ultra® Test Strip package. (Sold separately.)

WARNING: Keep the test strip vial away from children; the cap is a choking hazard. Also, the cap may contain drying agents that are harmful if inhaled or swallowed and may cause skin or eye irritation.

BEFORE TESTING

Checking the Display

Each time you turn on the OneTouch® Ultra® Meter either by inserting a test strip or pressing the M button, all segments of the display will appear briefly. This tells you that the system is performing several self-checks to confirm that the meter is working properly. To check that all display segments are working, hold the C button down.



Coding the Meter

Code numbers are used to calibrate the OneTouch® Ultra® Meter with OneTouch® Ultra® Test Strips for accurate results. You must code the meter before using it for the first time and then every time you change to another vial of test strips. **Each time you test, check to be sure that the code number on the meter display matches the code number on the test strip vial.**

CAUTION: Failure to code the meter correctly will cause inaccurate test results.


STEP 1 Enter the Code Mode.

Insert a test strip to turn on the meter. Push it all the way in until it will go no further. Avoid bending the test strip. The display check will appear. Then the code number is displayed for three seconds. (The first time you use the meter, three dashes - - - will appear, meaning that there is no code stored in the memory.) If three dashes appear any other time, see page 48 of “Display Messages and Problem-Solving Guide.”

STEP 2 Match the Code Numbers.

Compare the code number on the meter display with the code number on the test strip vial. If the two code numbers match, you may begin testing. If they do not match, follow Step 3.

STEP 3 Code the Meter.

Press the C button to select the correct code. Each time you press and release the C button, the number will increase by one. To move more quickly, press and hold the C button.

After selecting the correct code number, it will flash for three seconds and then appear solid for three seconds.

Then the  symbol will appear with the unit of measure, indicating that the OneTouch® Ultra® System is ready for testing.

Checking the System

NOTE: Refer to additional information in the OneTouch® Ultra® Control Solution package.

OneTouch® Ultra® Control Solution is used to check that the meter and the test strips are working together as a system and that you are performing the test correctly.

OneTouch® Ultra® Control Solution contains a measured amount of glucose that reacts with OneTouch® Ultra® Test Strips. Compare your control solution test results with the range printed on the test strip vial label. It is very important that you do this simple check routinely to make sure you get an accurate result.

Before you use the OneTouch® Ultra® Meter to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the range, you are ready to test your blood.

Important Control Solution Test Information

- Use only OneTouch® Ultra® Control Solution.
- Check the expiration date on the control solution vial. Record the discard date (date opened plus three months) on the vial label. Do not use after expiration or discard date, whichever comes first.
- Control solution, meter, and test strips should be at room temperature (20–25°C / 68–77°F) before testing.
- Shake the vial, discard the first drop of control solution, and wipe off the tip to ensure a good sample and an accurate result.
- Store control solution tightly closed at temperatures below 30°C (86°F). Do not refrigerate.

CAUTION: The control solution range printed on the test strip vial is for OneTouch® Ultra® Control Solution only. **It is not a recommended range for your blood glucose level.**

When to do a control solution test:

- Once a week.
- When you open a new vial of test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- If your blood glucose test results are not consistent with how you feel.
- After dropping the meter.

How to do a control solution test:


STEP 1 Insert Test Strip.

Insert a test strip, contact bars end first and facing up, into the test port. Push it all the way in until it will go no further. The meter will turn on and the display check will appear briefly.

Then the code number will appear, followed by the  symbol and unit of measure. Check that the unit of measure is set correctly. Be sure the meter and test strip codes match. If they do not, code the meter correctly.



Press the C button to mark the test as a control solution test in the meter memory. **ctl** will appear on the display. If you decide not to do a control solution test, press the C button again to remove **ctl** from the display.

STEP 2 Apply Control Solution.

To ensure an accurate result:


- Shake the vial well
- Discard the first drop
- Wipe the dispenser tip

Hold the drop to the narrow channel in the top edge of the test strip. When the confirmation window is full, the meter will begin to count down from 5 to 1 second.

STEP 3 Result Appears in 5 Seconds.

Compare the control solution test result with the range printed on the test strip vial. The result should fall within this range.

NOTE: Mark all control solution tests with **ctl** to distinguish them from blood glucose tests in the meter memory. Marked control solution tests will not be included in your averages.

Comparing control solution results

If test results fall outside the range printed on the test strip vial, repeat the test. Out-of-range results may be caused by one or more of the following:

- Error in performing the test.
- Failure to shake the control solution vial well.
- Expired or contaminated control solution.
- The meter, test strips, or control solution are too warm or too cold.
- Failure to discard the first drop of control solution and wipe the dispenser tip clean.
- Improper meter coding.
- Test strip deterioration.
- Meter malfunction.

CAUTION: If you continue to get control solution test results that fall outside of the range printed on the vial, the system may not be working properly. **Do not** use the meter. Call the OneTouch® Customer Care Line at **1 800 663-5521**.


TESTING YOUR BLOOD

Read this section and the test strip package insert carefully before testing. Make sure you have all items needed to test.

- Meter
- Test Strips
- Sampler
- Optional OneTouch® UltraClear™ Cap (for forearm sampling)
- Sterile Lancet

Preparing the OneTouch® UltraSoft™ Sampler


CAUTION: To reduce the chance of infection:


- Never share a lancet or the OneTouch® UltraSoft™ Sampler with anyone.
- Always use a new, sterile lancet. Lancets are for single use only.
- Keep the OneTouch® UltraSoft™ Sampler clean.

STEP 1 Insert a Lancet.


Turn the cap counterclockwise to remove it.


Insert the lancet into the lancet holder and push down firmly until it is fully seated. Do not twist the lancet. Twist the protective disk until it separates from the lancet. Replace the OneTouch® UltraSoft™ Cap. Turn it clockwise until it is snug.


Adjust the puncture depth setting if necessary. Twist the depth adjustment knob toward the smaller bumps for a shallower puncture or toward the larger bumps for a deeper puncture.

STEP 2 Cock the Sampler.

Slide the ejection/cocking control back until it clicks. If it does not click, the sampler may have been cocked when the lancet was inserted. The sampler is now ready for use.

STEP 3 Wash Your Hands and the Puncture Site.


Use warm, soapy water. Rinse and dry thoroughly.

Fingertip Blood Sampling

OneTouch® Ultra® System requires a very small blood drop to perform a test. You may obtain it from a **fingertip** or **forearm**. (See page 26 for information on obtaining a blood sample from the forearm.) Choose a different puncture site each time you test. Repeated punctures in the same spot may cause soreness and calluses.


STEP 1 Position the Sampler.

Hold the OneTouch® UltraSoft™ Sampler **firmly** against the **side** of your finger. Press the release button.


STEP 2 Massage the Fingertip.

Massaging the fingertip gently will help you obtain a round drop of blood. Do not squeeze excessively on the puncture site. The blood sample must be at least one microlitre (1 μ L) in volume (● actual size) or you may get an **Er 5** message (see page 55) or an inaccurate test result. Do not smear the blood sample. Proceed with your blood glucose test.


Forearm Blood Sampling

The forearm has fewer nerve endings than the fingertip so you may find that obtaining a blood sample from the forearm may be much less painful than using the fingertip. The technique for forearm sampling is different from fingertip sampling. Also there are differences between forearm samples and fingertip samples that you should understand. Please read the important information on page 28 carefully.

STEP 1 Install the OneTouch® UltraClear™ Cap.

To aid in obtaining a blood sample from the forearm, replace the regular sampler cap with the UltraClear™ Cap. If necessary, set the sampler for a deeper puncture.


STEP 2 Choose the Puncture Site.

Select a soft, fleshy area on your forearm that is clean and dry, away from bone, and free of visible veins and hair.


STEP 3 Massage the Area.

To increase blood flow to the puncture site, massage the area gently. For individuals who experience difficulty in getting sufficient blood for a test, rubbing the area more vigorously or applying heat briefly may be helpful.

STEP 4 Position the Sampler.

Press and hold the sampler against the forearm for a few seconds. Press the release button.

STEP 5 Allow Blood Drop to Form.

Continue holding the sampler against the skin for a few seconds until the blood drop forms. Allow enough blood to form under the cap until you have a blood sample that is sufficient to fill the confirmation window of the test strip. If you must massage the area to obtain more blood, do not squeeze the site excessively.

Important Information About Forearm Testing

- Under certain conditions, blood glucose test results obtained using samples taken from your forearm may differ significantly from fingertip samples.
- The conditions in which these differences are more likely to occur are when your blood glucose is changing rapidly such as following a meal, an insulin dose, or associated with physical exercise.
- When blood glucose is changing rapidly, fingertip samples show these changes more quickly than forearm samples.
- When your blood glucose is falling, testing with a fingertip sample may identify a hypoglycemic (low blood sugar) level sooner than a test with an forearm sample.
- Use forearm samples only for testing prior to, or more than two hours after meals, insulin doses, or physical exercise.
- Testing performed within two hours after meals, insulin doses, or physical exercise, or whenever you feel that your glucose levels may be changing rapidly, should be done with a fingertip sample.
- You should also use fingertip testing whenever you have a concern about hypoglycemia (insulin reactions) such as when driving a car, particularly if you suffer from hypoglycemic unawareness (lack of symptoms to indicate an insulin reaction), as forearm testing may fail to detect hypoglycemia.

Important Information About Forearm Testing (continued)

What you should do:

- Use either forearm or fingertip samples for testing prior to, or more than two hours after, meals, insulin doses, or physical exercise.
- Routine testing before meals can be done either at the fingertip or the forearm.
- Consult your healthcare professional before you begin using the forearm for testing.


If bruising occurs, you may choose to lance a fingertip instead. If you are having difficulty obtaining blood from the forearm, call the OneTouch® Customer Care Line at **1 800 663-5521**.

Step-by-Step Test Procedure


STEP 1 Insert Test Strip.

Insert a test strip, contact bars end first and facing up, into the test port. Push it in until it will go no further. The meter will turn on and the display check will appear briefly. Then the code number will appear, followed by the  symbol with the unit of measure. Check that the unit of measure is set correctly. Be sure the meter and test strip codes match. If they do not, code the meter correctly.


STEP 2 Apply Sample.

Obtain a round drop of blood using the OneTouch® UltraSoft™ Adjustable Blood Sampler. **The blood sample must be at least 1 μL in volume** () to fill the confirmation window. When the  symbol appears on the display, **touch and hold the drop of blood to the narrow channel in the top edge of the test strip.**



DO NOT

- Apply sample to the front or back of the test strip.
- Push your finger against the test strip.
- Apply a smeared sample.


Hold the blood drop to the top edge of the test strip until the confirmation window is full before the meter begins to count down.


If the confirmation window does not fill completely before the meter begins to count down, do not add more blood to the test strip; discard the test strip and retest. If you have trouble filling the test strip, call the OneTouch® Customer Care Line at **1 800 663-5521** for assistance.

If the confirmation window is not full, you may get an **Er5** message or an inaccurate test result.

NOTE: If you do not apply a blood sample within two minutes after the  symbol appears, the meter will turn itself off. You must remove the test strip and insert it back into the meter to restart the test procedure.


STEP 3 Accurate Results in Just 5 Seconds

Your blood glucose test result will appear after the meter counts down from 5 to 1. Blood glucose test results are automatically stored in the meter memory. You may also record them in the logbook provided with your system. Turn the meter off by removing the test strip. Dispose of used test strip in a sealed container.


Used Lancet Disposal

CAUTION: It is important to discard the used lancet carefully.

STEP 1 Remove the Sampler Cap.

Twist the cap counterclockwise. (Optional: Replace the protective disk on the used lancet by placing it on a hard surface and pushing the lancet tip into the disk.)

STEP 2 Eject the Lancet.




Point the sampler down and away. Push forward on the ejection/cocking control and eject the lancet directly into a container for sharp objects. Return the ejection/cocking control to the middle position. Replace the cap.

Special Messages

The OneTouch® Ultra® Meter displays results between 1.1 to 33.3 mmol/L (20 and 600 mg/dL). If your blood glucose test result is lower than 1.1 mmol/L (20 mg/dL), **LO** will appear on the meter display. **This indicates severe hypoglycemia (low blood glucose). You should treat this condition according to the recommendations of your healthcare professional.**


If your blood glucose test result is above 33.3 mmol/L (600 mg/dL), **H I** will appear on the meter display. **This indicates severe hyperglycemia (high blood glucose). You should recheck your blood glucose level. If H I again appears on the display, call your doctor immediately.**

When your blood glucose test result is above 13.3 mmol/L (240 mg/dL) or reads **H I** “**ketones?**” will appear on the meter display. **This message does not mean that the system detected ketones but that testing with a ketone test strip may be advisable. Consult your healthcare professional about when and how to test for ketones.**


USING THE METER MEMORY

Your OneTouch® Ultra® Meter stores the 150 most recent blood glucose and control solution test results with date and time in its memory. It also provides you with 14- and 30-day averages of your blood glucose test results. Review the test results in memory with these easy steps.

STEP 1 Enter the Memory Mode.


You may enter the memory mode either with the meter turned off or immediately after completing a test. Press the M button. **mem** and the 14-day average will appear with the number of blood glucose tests performed in this period. After three seconds, the 30-day average will appear. (If you are using the meter for the first time, - - - indicates there are no test results in memory.)

These averages are calculated from the blood glucose results obtained during the last 14 and 30 calendar days. A **H I** result will be included in your averages as 33.3 mmol/L (600 mg/dL) and a **L 0** result as 1.1 mmol/L (20 mg/dL). Results marked as control solution tests will not be included in your averages.


STEP 2 Recall Test Results.

The most recent test result with date and time will appear. Press and release the C button and the next most recent test result will appear. The meter will recall up to your last 150 test results in order from the most recent to the oldest. When the memory is full, the oldest result is dropped and the newest is added. To move through the tests more quickly, press and hold the C button. **ctl** will appear with control solution tests.

STEP 3 Exit the Memory Mode.

Press the M button to turn off the meter.


Downloading Test Results to a Personal Computer


Transferring test results with date and time from the meter memory to your personal computer for home viewing requires OneTouch™ Diabetes Management Software from LifeScan and a LifeScan Interface Cable both of which are not included with your OneTouch® Ultra® System kit. Visit www.OneTouch.ca to download the OneTouch™ Diabetes Management Software. The LifeScan Interface Cable can be purchased separately by calling the OneTouch® Customer Care Line at 1 800 663-5521.

STEP 1 Install OneTouch™ Software.

Follow the instructions provided in the OneTouch™ Software User's Manual.

STEP 2 Connect the LifeScan Interface Cable.

Connect the LifeScan Interface Cable to a serial port on the back of your computer and (with the meter turned off) to the meter data port.  will appear on the meter display when the first command is received, indicating that the meter is in the communication mode.


NOTE: While in the communication mode, you will be unable to perform a blood glucose test. If the meter is in either memory or setting mode, it will not respond to computer commands.

STEP 3 Transfer Data.

Follow the instructions in the OneTouch™ Software User's Manual to initiate the data transfer command. If the command is not received within two minutes, the meter will turn itself off. Pressing the M button will also turn the meter off.

To learn more about OneTouch™ Diabetes Management Software or to obtain a LifeScan Interface Cable separately, contact the OneTouch® Customer Care Line at **1 800 663-5521**.

COMPARING METER AND LABORATORY RESULTS

OneTouch® Ultra® Meter test results and laboratory test results are both expressed in plasma-equivalent units. However, your meter result may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. Your OneTouch® Ultra® Meter glucose value should agree with a laboratory measurement within $\pm 20\%$ most of the time under normal conditions. A result within that range is considered accurate when testing with the OneTouch® Ultra® Meter. However, results can differ by more than $\pm 20\%$ in some situations. See the OneTouch® Ultra® Test Strip package insert for typical accuracy and precision data and for important information on Limitations of Procedure.

To make an accurate comparison between meter and laboratory results, follow these guidelines:

Before you go to the lab:

- Perform a control solution test to make sure the meter is working properly.
- It is best to fast for at least eight hours before doing comparison tests.
- Take your meter with you to the lab.

While at the lab:

- Make sure that the meter test and the lab test are performed within 15 minutes of each other.
- Wash your hands before obtaining a blood sample.
- Use only fresh capillary blood obtained from the fingertip.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication, or experienced stress.¹ In addition, if you have eaten recently, the blood glucose level from a fingerstick can be up to 3.9 mmol/L (70 mg/dL) higher than blood drawn from a vein (venous sample) used for a lab test.² Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a meter result to be different from a laboratory result.

References

1. Surwit, R.S., and Feinglos, M.N.: *Diabetes Forecast* (1988), April, 49–51.
2. Sacks, D.B.: “Carbohydrates.” Burtis, C.A., and Ashwood, E.R. (ed.), *Tietz Textbook of Clinical Chemistry*. Philadelphia: W.B. Saunders Company (1994), 959.

Setting the Meter Time, Date, and the Unit of Measure

The OneTouch® Ultra® Meter comes with the time and date preset. Check that the time, date, and the unit of measure are set correctly. If you need to change the time, date, or unit of measure; follow the instructions for changing the settings or call the OneTouch® Customer Care Line at 1 800 663-5521. You may also need to enter the setting mode if you remove the battery.

To set the time, date, or unit of measure, you must enter the setting mode. Start with the meter off. Then press and hold the M button for three seconds. The meter is now in the setting mode.

STEP 1 Set the Hour.

With the **hour** setting flashing, press and release the C button to advance one **hour**. To move faster, hold the C button down. With the correct hour on the display, press the M button and the minutes setting will start flashing.

Before setting the date, you must have first set the time.

STEP 2 Set the Minutes.

Press and release the C button to advance one **minute**. To move faster, hold the C button down. With the correct minute on the display, press the M button and the AM/PM (12-hour) format setting will start flashing.

STEP 3 Set the Time Format.


The meter can display the time in either an AM/PM (12-hour) format or a 24:00 (24-hour) format. The AM/PM format is preset in the meter. If you wish to display the time in the 24:00 format, press the C button. With the preferred time format on the display, press the M button and the year setting will start flashing.

STEP 4 Set the Year.

Press and release the C button to advance one **year**. To move faster, hold the C button down. With the correct year on the display, press the M button and the date will appear on the display with the month segment flashing.

STEP 5 Set the Month.

Press and release the C button until the correct **month** appears. To move faster, hold the C button down. With the correct month on the display, press the M button and the day segment will start flashing.

STEP 6 Set the Day.

Press and release the C button until the correct **day** appears. To move faster, hold the C button down. With the correct day on the display, press the M button and the current unit of measure will start flashing.

Note: *The 14- and 30-day averages in the meter memory are calculated from results obtained during the 14 and 30 calendar days preceding the current date and time settings. When the date and time are changed, the averages may change.*

Note: *You must move through the hour, minutes, time format, year, month and day, to turn off the meter and exit the setting mode.*

CARING FOR YOUR OneTouch® Ultra® SYSTEM

Meter

Your OneTouch® Ultra® Meter does not require special maintenance or cleaning. Take care to avoid getting dirt, dust, blood, control solution, or liquids inside the meter through the test port or data port. LifeScan recommends that you store the meter in its carrying case after each use.

A cloth dampened with water and mild detergent can be used to wipe down the outside of the meter. Your OneTouch® Ultra® Meter is a precision instrument. Please handle it with care.

OneTouch® UltraSoft™ Sampler



Clean the sampler and caps with soap and water. To disinfect the sampler, prepare a disinfectant solution of one part household bleach to 10 parts water. Dampen a cloth with this solution and wipe the sampler thoroughly. **Soak only the cap or caps** for at least 30 minutes in the disinfectant solution. Do not soak the sampler in liquid. Rinse the sampler and cap with water and dry thoroughly.

Battery


Your OneTouch® Ultra® Meter comes with a pre-installed 3.0 V (#2032 or equivalent) lithium battery. The battery will provide you with enough power to perform about 1,000 tests. When replacing the battery, use only a 3.0 V (#2032 or equivalent) lithium battery.

The meter will alert you when the power is getting low by displaying two different messages:

1. The  symbol appears on the display with the unit of measure when the meter is turned on and all other display messages are functional. From the time the  symbol first appears, there is enough power left for about 50 tests. The test results will be accurate, but it is time to change the battery.
2. When the  symbol appears on the display by itself, it means that the battery will not provide enough power for a test. You must change the battery.


To replace the battery, make sure that the meter is turned off. Turn the meter over and locate the battery compartment.


STEP 1 Open the battery compartment.

Push up on the tab and pull the battery compartment cover toward you.

STEP 2 Remove the old battery.

Pull upward on the ribbon. Insert one 3.0 V (#2032 or equivalent) lithium battery, making sure the positive “+” side of the battery is facing up.

STEP 3 Place the cover.

Place the cover over the battery compartment. Push down until you hear the tab click into place.


Note:

- Replacing the battery **does not** affect previous test results stored in memory. However, the time and date settings may need to be updated.
- After replacing the battery, turn the meter on by inserting a test strip or pressing the M button. If the meter is in setting mode, confirm that the time, time format and date are set correctly. If they are not, use the M and C buttons to reset the meter before testing. See “Setting the Meter,” pages 40–43.

DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE

Following is a summary of all display messages. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message or a symbol. In the event of a problem, refer to the information under ACTION. Assistance is available from the OneTouch® Customer Care Line.

MESSAGE


WHAT IT MEANS

Display check. Appears whenever the meter is turned on. To verify that all segments are present on the display, press and hold the C button. Compare the meter display with page 12.




ACTION


If segments are missing, contact the OneTouch® Customer Care Line for further action. Missing segments can lead to misinterpretation of test results.

- 1) This message appears immediately after the display check if your meter has not been coded.
- 2) If this message appears at any time after the first time meter was coded, the code number has been lost. The test results stored in the meter memory may be out of order.

- 1) See “Coding the Meter,” pages 12-14.
- 2) Do not rely on the 14- and 30-day averages. It is important to call the OneTouch® Customer Care Line.

<u>MESSAGE</u>	<u>WHAT IT MEANS</u>	<u>ACTION</u>
 <p>The image shows a square frame containing the word "code" in the top left corner and the number "12" in a large digital font in the center.</p>	<p>This is the code number stored in the meter.</p>	<p>Make sure that this code number matches the code number on the vial of test strips you are using.</p>
 <p>The image shows a square frame containing a blood drop icon above a horizontal line, and the unit "mmol/L" below it.</p>	<p>The system is ready to accept a blood sample.</p>	<p>You may now apply the blood sample.</p>
 <p>The image shows a square frame containing the label "ctl" in the top left corner, a control drop icon above a horizontal line, and the unit "mmol/L" below it.</p>	<p>The system is ready to accept a control solution test sample.</p>	<p>You may now apply the control solution.</p>
 <p>The image shows a square frame containing the number "5" in a large digital font.</p>	<p>Five-second countdown. The meter is calculating the result. After the countdown, the meter will display the test result.</p>	<p>None required.</p>

MESSAGE**WHAT IT MEANS**

A blood glucose test result in mmol/L.

ACTION


None required.

A blood glucose test result with a suggestion to check your ketone levels.

Follow the instructions of your healthcare professional regarding ketone testing.

MESSAGE

WHAT IT MEANS

You may have a very high blood glucose level, exceeding 33.3 mmol/L (600 mg/dL).


ACTION

You should recheck your glucose level and if HI call your doctor immediately.

You may have a very low blood glucose level, lower than 1.1 mmol/L (20 mg/dL).


Treat this condition according to your healthcare professional's recommendations.

A blood glucose test result stored in the memory.

None required.

MESSAGE


(Example)

WHAT IT MEANS

A control solution test result stored in the memory.


ACTION

None required.

There is no test result stored in this place in the memory.


None required.

Average of the last 14 days 64_n means that 64 tests were performed within the period). The 30-day average appears with a 30 .

None required.

MESSAGE


WHAT IT MEANS

Error message that indicates that there is a problem with the meter.

ACTION

Do not use the meter. Contact the OneTouch® Customer Care Line at **1 800 663-5521** for a replacement meter.

Error message could be caused either by a used test strip or a problem with the meter.

Repeat the test with a new test strip. Refer to page 10 for test strip information. If the error message persists, contact the OneTouch® Customer Care Line at **1 800 663-5521**.

Error message that indicates that the blood or control solution sample was applied before the  symbol appeared on the display.


Repeat the test with a new test strip. Apply blood or control solution sample only after the  symbol appears on the display. Refer to pages 30-31 for information on sample application.

MESSAGE

WHAT IT MEANS

ACTION

Error message that indicates one of the following conditions may be present:

1) You may have high glucose and have tested in an environment near the low end of the system's operating temperature range(6–44°C/ 43–111°F).

If you tested in a cool environment, repeat the test in a warmer environment with a new test strip.*

2) There may be a problem with the test strip. For example, it may have been damaged or moved during testing.


If you were testing in a normal or warm environment, repeat the test with a new test strip. (See pages 10-11 for test strip information.)*

3) The sample was improperly applied.

If you applied the blood incorrectly, review pages 30-31 on blood application and repeat the test with a new test strip.*

*If the error message appears again, contact OneTouch® Customer Care Line at **1 800 663-5521.**

MESSAGE

WHAT IT MEANS

Error message that indicates that the meter has detected a problem with the test strip. Possible causes are test strip damage or an incompletely filled confirmation window.


ACTION

Repeat the test with a new test strip. Refer to pages 30-31 for information on sample application.

The  symbol appears on the display with the unit of measure. The power is getting low. You can complete about 50 more tests from the time this symbol first appears.

Test results will still be accurate, but replace the battery as soon as possible.

The  symbol appears on the display by itself. The power is too low to perform a test.

Replace the battery at once. The meter will not operate.

MESSAGE**WHAT IT MEANS**

Flasing time means the meter has lost power and entered the setting mode.

ACTION

You must reset the time, time format and date. See “Setting the Meter,” pages 40-43.

The meter has detected that the temperature is below the system operating range. Do not perform a test until the meter and test strips reach a temperature within the operating range of 6–44°C (43–111°F).

Repeat the test after the meter and test strips have reached a temperature within the operating range.

The meter has detected that the temperature is above the system operating range. Do not perform a test until the meter and test strips reach a temperature within the operating range of 6–44°C (43–111°F).

Repeat the test after the meter and test strips have reached a temperature within the operating range.

If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE

WHAT TO DO

Need battery

Replace the battery.

No battery or battery incorrectly installed

Check that the battery is correctly installed with the positive “+” side up.

Test strip inserted incorrectly or incompletely

Insert the test strip correctly with the contact bars end first and facing up. Push it in until it will go no further.

Defective meter or test strip

Call **1 800 663-5521**.

If the test does not start after applying the sample:

POSSIBLE CAUSE

WHAT TO DO

Insufficient blood sample

Repeat the test with a new test strip and a larger sample.

Defective test strip

Repeat the test with a new test strip.

Sample applied after automatic shutoff (two minutes after last user action)

Repeat the test with a new test strip; apply sample only when  appears on the display.

Defective meter

Call **1 800 663-5521**.

SPECIFICATIONS

Result Range:	1.1 to 33.3 mmol/L (20 to 600 mg/dL)
Calibration:	Plasma-equivalent
Sample:	Fresh capillary whole blood
Sample Size:	Minimum 1 microlitre
Test Time:	5 seconds
Assay Method:	Glucose oxidase biosensor
Power Source:	One replaceable 3.0 V (#2032 or equivalent) lithium battery
Battery Life:	1,000 tests or about one year at three tests per day
Glucose Units:	mmol/L
Memory:	150 blood glucose and control solution tests
Automatic Shutoff:	Two minutes after last user action
Size:	7.9 cm x 5.7 cm x 1.9 cm
Approximate Weight:	42.5 grams with battery
Operating Ranges:	Temperature 6–44°C/43–111°F Relative Humidity 10–90% Hematocrit 30–55%

GUARANTEE

The makers guarantee that the OneTouch® Ultra® Meter shall be free of defects in material and workmanship for a period of three years. This guarantee is valid from the date of purchase. The guarantee extends only to the original purchaser and is not transferable.

The OneTouch® Ultra® Blood Glucose Meter complies with the following standards:

CAN/CSA C 22.2 No. 601.1 – M90

“Medical Electrical Equipment –

Part 1: General Requirements for Safety” (Complies with UL 2601-1)

EN 60601-1-2 (1993) “Medical electrical equipment: Part 1, Section 1.2 Collateral standard:”
EMC requirements & tests.”

EN 50082-1 (1998) “Electromagnetic Compatibility, Generic Immunity Standard, Part 1, Residential, Commercial & Light Industry.”

- Internally powered equipment.
- No patient-applied parts.
- This unit is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

INDEX

Averages, 14- and 30-Day.....	34
Battery.....	45
Blood Sampling, Forearm.....	26
Blood Sampling, Fingertip.....	25
C Button.....	7
Checking the Meter Display.....	12
Checking the System.....	15
Coding the Meter.....	12
Control Solution, Checking with.....	15
Control Solution, OneTouch® Ultra®.....	15
Control Solution Results, Comparing.....	20
Control Solution Test Information, Important.....	16
Data Downloading.....	36
Data Port.....	36
Display Messages/Problem-Solving Guide.....	48
Guarantee.....	59
Laboratory Results, Comparing Meter Results with.....	38
Lancet.....	23
Lancet Disposal.....	32
LifeScan Interface Cable.....	36
M Button.....	6
Memory, Using.....	34
Messages, Special.....	33

Meter Display	8
Meter Display, Checking.....	12
Meter Maintenance	44
Meter Specifications.....	58
OneTouch™ Software	36
OneTouch® Ultra® Blood Glucose Meter	6
OneTouch® Ultra® Blood Glucose Monitoring System	5
OneTouch® Ultra® Control Solution	15
OneTouch® Ultra® Test Strips	10
- Code.....	12
- Expiration.....	11
OneTouch® UltraSoft™ Sampler, Cleaning	44
OneTouch® UltraSoft™ Sampler, Preparing	22
OneTouch® UltraSoft™ Sampler, Using.....	23
Plasma Calibration.....	iv
Problems, Solving	48
Serial Port	36
Setting the Time and Date	40
Step-by-Step Test Procedure	30
System, Checking the	15
Test Port	7
Test Procedure, Step-by-Step	30
Test Strip Information, Important	11
Test Strips, OneTouch® Ultra®.....	10

Times of Day	Glucose Ranges for People Without Diabetes, mmol/L (mg/dL)	Your Target Ranges, mmol/L (mg/dL)
Before breakfast	3.9–5.8 (70–105)	_____
Before lunch or dinner	3.9–6.1 (70–110)	_____
1 hour after meals	Less than 8.9 (160)	_____
2 hours after meals	Less than 6.7 (120)	_____
Between 2 and 4 AM	Greater than 3.9 (70)	_____

Source: Krall, L.P., and Beaser, R.S.: *Joslin Diabetes Manual*. Philadelphia: Lea and Febiger (1989), 138.

Meter Serial No. _____

Important phone numbers:

OneTouch® Customer Care Line (Monday–Friday 9AM–8PM Eastern Time, 6AM–5PM Pacific Time): 1 800 663-5521

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©LifeScan, Inc. 2003 Milpitas, CA 95035 02/2004