
Ear Thermometer

GRACO

OWNER'S MANUAL

KEEP INSTRUCTIONS FOR FUTURE REFERENCE.

©2008 Graco PD104815A 9/08

WARNING

Failure to follow these warnings and the assembly instructions could result in serious injury or death.

- This product not to be used as a toy.
- **WARNING:** This product is not a toy.
- **DO NOT** immerse any part of the Ear Thermomete, in water.
Clean only with a dry cloth.
- **DO NOT** place unit near water and moisture. Do not use the Ear Thermometer outdoors. Do not use the Ear Thermometer near possible wet areas, such as a bathtub, shower, wash basin, sink, laundry tub, swimming pool, wet basement, etc.
- **Keep Ear Thermometer away from heat sources**, such as fireplaces, radiators, stoves, and TVs or other appliances. Heat can damage the case or electrical parts.

• **DO NOT** use the Ear Thermometer Unit if:

- The unit has been exposed to liquids.
- Unit has been dropped or is damaged.

• **DO NOT OPEN THE EAR THERMOMETER.**

No user-serviceable parts inside. Risk of electrical shock, fire or death.

- Do not drop the instrument! Protect it from severe impact and shock.

FOR SAFE BATTERY USE:

Keep the batteries out of children's reach. Any battery may leak battery acid if mixed with a different battery type, if inserted incorrectly (put in backwards) or if all batteries are not replaced or recharged at the same time. Do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (nickel-cadmium) batteries.

Any battery may leak battery acid or explode if disposed of in fire or an attempt is made to charge a battery not intended to be recharged. Never recharge a cell of one type in a recharger made for a cell of another type.

Discard leaky batteries immediately. Leaky batteries can cause skin burns or other personal injury. When discarding batteries, be sure to dispose of them in the proper manner, according to your state and local regulations.

Remove batteries before storing product for a prolonged period of time. Batteries left in the unit may leak and cause damage.

Recommended type batteries—disposables, size 3-volt CR 2032 Lithium battery. NEVER mix battery types. Change batteries when product ceases to operate satisfactorily.

Always remove batteries if the product is not going to be used for a month or more. Batteries left in the unit may leak and cause damage. Recommended type batteries—disposables, size 3-volt CR 2032 Lithium battery. NEVER mix battery types. Change batteries when product ceases to operate satisfactorily.

WARNING

Modifications not authorized by the manufacturer may void users authority to operate this device.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:



- Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.


This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Overview

Accessories (On certain models)



Symbol Key

	Degrees of temperature
	Temperature is being scanned
	Smiling face indicating temperature lower than 100.4°F (38°C)
	Frowning face indicating temperature of 100.4°F (38°C) or higher
°F	Fahrenheit
°C	Celsius
	Low battery warning
	Replace battery
	Memory symbol
	Temperature is below 93.2°F (34°C)
	Temperature is above 109.4°F (43°C)
	Ambient temperature is outside thermometer's range 60.8°-104°F (16°-40°C)



Battery Installation

NOTE: The word “Lo” by itself does not mean that the battery is low and needs to be charged. It means that the temperature being measured is out of the thermometer’s range (see Fig. 8).

Replace the battery when the low battery symbol appears (see Fig. 1). If only the symbol appears, you may still use the thermometer but you will need to replace the battery soon. If the symbol appears with the word “Lo” (see Fig. 2), you must replace the battery to use the thermometer.


1. Slide the battery compartment lid off.
2. Carefully push back on the battery release using a non-metal tool. Battery will pop up (see Fig. 4).
3. Remove old battery from compartment and discard immediately.
4. Insert one 3-volt CR 2032 lithium battery with the positive side (+) facing up.
5. Replace the battery compartment lid by sliding it forward.


Instructions for use

Please do not use this ear thermometer when the ear canal is wet.


1. Press and release the Power/Scan button. The system will start its self-testing cycle (see Fig. 5).


2. Once the self-testing is completed (approximately 1 second), the thermometer will beep and the display will show the last temperature recorded. This indicates the thermometer is ready to take a temperature (see Fig. 6).


3. Press the Power/Scan button again. The symbol  will appear beneath the last reading (see Fig. 7). The word "Lo" may also appear




indicating the temperature is lower than 93.2°F (34°C) (see Fig. 8). NOTE: The word "Lo" by itself does not mean that the battery is low and needs to be charged unless the low battery icon appears with it (see Fig. 2). It means that the temperature being measured is out of the thermometer's range.





4. Insert the probe into the ear canal. The probe must be inserted well into the ear to allow the probe tip to continuously scan the ear canal. This will ensure an accurate temperature measurement (see. Fig. 9).

- After one second you will hear a beep indicating the reading is complete. Remove the thermometer from the ear to get the temperature reading.
- A smiling face will appear indicating a temperature reading lower than 100.4°F (38°C). A frowning face accompanied by one long and two short beeps will indicate a temperature of 100.4°F (38°C) or higher (see Figs. 10 & 11).
- You do not need to manually shut off the thermometer. It will automatically shut off after one minute of non-use.
- To change the reading from Celsius to Fahrenheit, press the Power/Scan button without releasing it, then press the blue M button on the underside of the thermometer. First release the M button and then the Power/Scan button.
- The Ear Thermometer stores and recalls the last 10 temperature readings. To recall the prior readings simply press the blue M button on the underside of the thermometer. The memory symbol, the number of the reading, and the temperature reading will appear on the display (see Fig. 12).

10. If you want to immediately take another reading, please wait at least 10 seconds.
11. The temperature may vary between both ears so try to take it from the same ear each time.
12. The thermometer has a temperature range of 93.2°-109.4°F (34°-43°C). If the temperature being measured is below the range, the word “Lo” will appear. The word “Hi” will appear when the temperature measured is above 109.4°F (43°C) (see Figs. 13 &14)
13. If the ambient temperature is outside the thermometer’s range (60.8°-104°F/16°-40°C), the word “Err” will appear (see Fig. 15).

Care and maintenance ● ● ● ● ● ●

1. The ear thermometer is a delicate electronic device. Although the probe can be cleaned with isopropyl alcohol, the body of the thermometer should never be submerged in water or any other liquid.
2. Clean the probe with a cotton swab moistened with isopropyl alcohol after every measurement and before each use. Clean the body of the thermometer with a soft dry cloth. Never use an abrasive.
3. Avoid touching the probe with your fingers.
4. Do not expose the thermometer to direct sunlight.
5. Keep out of the reach of children when not in use. The optimal storage temperature range is 4°-122°F (-20°-50°C).
6. Use the ear thermometer to take temperature from the ear canal only. Do not use it to take rectal, oral, or underarm temperatures.
7. Do not disassemble the thermometer.
8. Only use a 3-volt CR-2032 lithium battery.

Troubleshooting

1. Consistent low temperature readings.
 - a. The probe is not positioned properly in the ear canal. The tip of the probe must be snug and fully seated against the opening of the ear canal. Failure to properly position the probe may lead to a low temperature reading.
 - b. The probe lens is dirty. Clean the lens with a soft cloth moistened with water or isopropyl alcohol.
2. The word "Lo" or "Hi" appears. Temperature is out of the range of thermometer.
 - a. The probe is not positioned properly in the ear canal. The tip of the probe must be snug and fully seated against the opening of the ear canal. Please properly position the ear thermometer in the ear canal and then repeat the measurement.
3. Thermometer will not take a reading.
 - a. Check to see if the low battery icon appears with the word "Lo". If it does, the battery needs to be replaced. See Battery Installation section of the manual.
 - b. Check to see if the low battery icon appears by itself. If it does, the battery may need to be replaced. See the Battery Installation section of the manual.

Applied standards and approvals ● ● ●

This product conforms to the provisions of the EC directive MDD(93/42/EEC).

The following standards apply to design and/or manufacturing of the products:

- ASTM E1965-98: Standard Specification for infrared thermometers for intermittent determination of patient temperature.
- IEC/EN980:1997: Graphical symbol for use in the labeling of medical devices.
- IEC/EN60601-1: Medical electrical equipment – Part 1: General requirement for safety.
- IEC/EN60601-1-2: Medical electrical equipment – Part 2: Collateral standard: Electromagnetic compatibility – Requirements and tests.
- ISO 14971 Risk Management

Type BF

Replacement Parts Warranty Information (USA)

To purchase parts or accessories or for warranty information in the United States, please contact us at the following:

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1-800-667-8184

Montreal (514-344-3533)

Fax: 514-344-9296


Product Registration (USA)

To register your Graco product from within the U.S.A. visit us online at

www.gracobaby.com/productregistration.

We currently do not accept product registrations from those living outside the United States of America.

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